



DoseWise: With every new system, Philips looks at how we can incorporate better shielding and improve X-ray exposure to further reduce dose. The Philips unique beam filters for example is a feature that reduces dose by 40% over conventional filters. And with the pulsed fluoro modes, dose can be reduced by up to 75%.



It's about reducing size...



Revolutionize the way you work with BV Endura. Optimally designed, the MobileView Station not only looks better, it performs better. Thanks to its compact size, this system can be easily moved around even in the tight confines of the OR. And its small footprint ensures that it can be brought as close to the operating table as is desired. Its unique design also enables easy cleaning of the equipment.

With the BV Endura and the intelligent viewing concept, the operator can enter patient demographics manually or retrieve a worklist via the hospital network. All this while the MobileView Station is already positioned at the surgical table. At the start of the procedure, the operator simply needs to turn the monitor towards the physician. Later, the monitors can be turned back to the operator to process the images and send them to PACS. To ensure high quality on-screen images, the MobileView Station is equipped with 18" LCD monitors. And to provide maximum flexibility, the height of the monitors can be adjusted to the surgeon's convenience – whether sitting or standing.

Height movement of the foldable monitors for optimal working and transport height.



and adding space.



The BV Endura is designed to give you more space in the OR by taking up less of it.

The compact MobileView Station is easy to transport, maneuver and makes more room for the surgeon and operator. In short, the BV Endura is not only designed to give you more space, it's designed to give you more freedom.





It's about optimizing workflow...



Increase your efficiency in many ways with the BV Endura. The BV Endura can be quickly moved into place, thanks to its ultra compact foot and C-arc handles. With the single user concept and Anatomically Programmed Fluoroscopy (APF) parameters, you can easily control every step of your procedures.



The BV Endura has a fully integrated DICOM connectivity solution (full RIS/HIS compatibility, export of examination onto the hospital network) and is equipped with a state-of-the-art ViewForum Surgical Workstation. With this workstation you can import multi-modality images directly into your workstation and view them instantly on the right hand monitor that acts as a digital interactive light-box.

An intuitive user-interface with buttons exactly where you need and expect them to be, makes working with the BV Endura easier and faster.

With a touch screen on the left monitor, patient administration and post-processing on acquired images is literally at your fingertips.

All extras like the paper/ transparency printer and Medical DVD recorder can of course be easily accommodated into your Mobile View Station.



Feel safe with Philips

Philips has over 50 years of experience in developing and servicing mobile C-arm systems. This proven history and experience is our solid basis to continue providing you with unparalleled customer service for the new BV Endura.

When you invest in the BV Endura, you invest in the future. Philips' ongoing commitments to develop future-safe technology means that you can keep your BV Endura system up-to-date throughout its lifecycle, embracing emerging application demands, and keeping up with advances in networking and PACS.

Through service excellence, flexible solutions and effective relationships, Philips' CUSTOMerCARE service programs give you all the support you need to ensure your BV Endura system always operates at peak performance – today and in the future.

and efficiency.



**Philips Medical Systems is part
of Royal Philips Electronics**

Interested?

Would you like to know more about our imaginative products? Please do not hesitate to contact us. We would be glad to hear from you.

On the web

www.medical.philips.com

Via email

medical@philips.com

By fax

+31 40 27 64 887

By postal service

Philips Medical Systems
Global Information Center
P.O. Box 1168
5602 BD Eindhoven
The Netherlands

BV Endura and CUSTOMerCARE are trademarks of Koninklijke Philips Electronics N.V. Certain features described in this brochure are optional or may not be available in some countries. Please contact your local Philips representative for further details.

Asia

Tel: +852 2821 5888

Europe, Middle East, Africa

Tel: +31 40 27 87246

Latin America

Tel: +55 954 628 1000

North America

Tel: +1-800-285-5585

©2005 Koninklijke Philips Electronics N.V.

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.

4522 962 08151/000 * OCT 2005

